

Jazzy[®]
SELECT 6



Colours Available:



Added Stability For Top Performance!

Sporty, Sophisticated and Unique; The Jazzy Select[®] 6 offers superior suspension for powerful performance and stability. Experience a tight turning radius and a depth-adjustable seat.

Features



Optional 4" Seat Lift Available



Intuitive Controller is Easy to Operate



Flip Out Foot Board



Active-Trac® Suspension

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length ³	8985 mm
Overall Width ³	696 mm (24.5")
Ground Clearance ³	50 mm (2.75")
Suspension	ATX Suspension (Active-Trac With Extra Stability)
Colour(s)	Blue & Red
Power	
Turning Radius ³	556 mm (20.5")
Max Speed ¹	4 mph
Drivetrain	Two Motor, Mid-Wheel-6 Drive
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements ³	2 x 12v, Deep Cycle NF-22 Batteries
Battery Weight ⁴	11.11 kg
Battery Charger	8 Amp Off-Board
Total Weight Without Batteries ^{3, 4}	45.36 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.