

Added Stability For Top Performance!

Sporty, Sophisticated and Unique; The Jazzy Select[®] 6 offers superior suspension for powerful performance and stability. Experience a tight turning radius and a depth-adjustable seat.





Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1067 mm
Overall Width	597 mm (23.8″)
Ground Clearance	19 mm (0.75″)
Suspension	ATX Suspension (Active-Trac With Extra Stability)
Colour(s)	Blue & Red
Power	
Turning Radius	546 mm (21.5″)
Max Speed	4 mph (6.4 km/h)
Drivetrain	Two Motor, Mid-Wheel-6 Drive
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements	2 x 12 Volt, Deep Cycle
Battery Weight	11.11 kg
Battery Charger	3.5 Amp Off-Board
Total Weight Without Batteries	45.3 kg (100 lbs)

Features



Optional 4" Seat Lift Available



Intuitive Controller is Easy to Operate



Flip Out Foot Board





www.pride-mobility.co.uk