

Jazzy[®] SELECT



Colours Available:



Highly Manoeuvrable Loaded With Performance And Convenience

The highly manoeuvrable Jazzy[®] Select delivers the advantages of in-line motor technology for enhanced efficiency, torque, range, and performance. Plus, it's loaded with performance and convenience features, like standard Active-Trac[®] Suspension, that make it easy to enjoy.

Features



Flip Out Foot Board



Intuitive Controller is Easy to Operate



Active-Trac® Suspension

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Max Speed	4 mph
Overall Length ⁴	90.80 cm
Overall Width ⁴	60.32 cm
Drive Tyres	25.4 cm (10")
Ground Clearance	4.76 cm (1.875") at Midframe
Suspension	Active-Trac®
Total weight (without batteries)	40.82 kg
Colour(s)	Red & Blue
Power	
Turning Radius ⁴	50.8 cm
Drivetrain	Two - Motor, Mid - Wheel Drive
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements ³	2 x 12v, Deep Cycle
Battery Weight ²	11.11 kg
Battery Charger	5 - Amp Off-Board

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.