

MAXIMA



Colours Available:



Ultra Heavy Duty With Maximum Performance

Heavy duty and robust, the Maxima scooter maximises your potential! Drive at speeds up to 5.8 mph and experience a fully adjustable tiller. Enjoy the convenience of the front basket, perfect for storing your personal items.

Specifications		
Scooter	4-Wheel	3-Wheel
Weight Capacity	36 st (227 kg)	
Overall Length	1206.5 mm (47.5")	
Overall Width	685.8 mm	
Drive Tyres	226.7 mm x 91.44 mm Solid 273.05 mm x 91.44 mm Solid	254 mm x 91.44 mm Solid 273.05 mm x 91.44 mm Solid
Ground Clearance	63.5 mm (2.5")	57.15 mm (2.25")
Colour(s)	Ocean Blue (Matte), Garnet Red (Matte)	
Mobility		
Turning Radius	1245 mm (49.01")	1346 mm (52.9")
Max Speed	Up to 5.8 mph	
Seat-to-Ground Height With Power Elevating Seat	561.34 mm - 612.14 mm	58.93 mm - 725.68 mm with Power Elevating Seat
Seat-to-Deck Height With Power Elevating Seat	435 mm - 485.9 mm	440.94 mm - 600 mm
Drivetrain	24v Single Motor	
Dual Breaking System	Regenerative & Electromechanical	
Battery		
Battery Requirements	Size: (2) Gel, 12v, 58 Ah, NF-22 Weight:16 55 kg. Each	Size: (2) Gel, 12v, 58 Ah, NF-22 Weight:16 55 kg. Each
Battery Charger	5A Off-Board	5A Off-Board
Total Weight Without Batteries	83.23 kg	89.26 kg


Features



Easy Control Tiller



Available in a 3 or 4 Wheeler



Comfortable Medium - Back Seat