

MAXIMA



Colours Available:



Ultra Heavy Duty With Maximum Performance

Heavy duty and robust, the Maxima scooter maximises your potential! Drive at speeds up to 5.8 mph and experience a fully adjustable tiller. Enjoy the convenience of the front basket, perfect for storing your personal items.


Features



Available in a 3 or 4 Wheeler



Easy Control Tiller



Comfortable Medium - Back Seat

Specifications		
Scooter	3 Wheel	4 Wheel
Weight Capacity	35st 10lbs 226.77 kg	
Overall Length ²	1206.5 mm	
Overall Width ²	685.8 mm	
Drive Tyres Back & Front ²	266.7 mm x 91.44 mm Solid 273.05 mm x 91.44 mm Solid	254 mm x 91.44 mm Solid 273.05 mm x 91.44 mm Solid
Ground Clearance ²	63.5 mm at Motor	57.15 mm at Motor
Colour(s)	Ocean Blue (Matte), Garnet Red (Matte)	
Mobility		
Turning Radius ²	1245 mm	1346 mm
Max Speed ¹	Up to 5.8 mph	
Seat-to-Ground Height ² with Power Elevating Seat	561.34 mm – 612.14 mm	58.93 mm – 725.68 mm with Power Elevating Seat
Seat-to-Deck Height ² with Power Elevating Seat	435 mm – 485.9 mm	440.94 mm – 600 mm
Drivetrain	24v Single Motor	
Dual Breaking System	Regenerative & Electromechanical	
Battery		
Battery Requirements ^{3, 5}	Size: (2) Gel, 12v, 58 Ah, NF-22 Weight: 16 kg. Each	
Battery Charger	5 - Amp Off-Board	
Total Weight Without Batteries ⁶	83.23 kg	89.26 kg

NOTE:

- Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- Includes standard seat weight.
- Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- Battery weight may vary based on battery manufacturer.
- Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

