

GO GO
ENDURANCE Li



Colours Available:



Travelling Has Never Been Easier!

Enjoy a Smooth, Comfortable Ride

With the Go-Go® Endurance Li, travelling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8 Ah airline-compliant lithium battery pack**, the Go-Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need when you arrive at your destination. Enjoy a smooth, comfortable ride with Comfort-Trac Suspension.

Features



Dual Charging Ports Can Charge Battery Pack On/Off-board.



Delta Tiller With Ergonomic Wraparound Handles



Standard 8 Ah Airline-Compliant Lithium Battery

Specifications	
Scooter	
Weight Capacity	23 st (147 kg)
Overall Length ⁴	1066 mm (42")
Overall Width ⁴	508 mm (20")
Drive Tyres ⁴	Front: 51 x 203 mm (2" x 8") - Solid Rear: 63 x 203 mm (2.5" x 8") - Solid
Ground Clearance ⁴	83 mm (3.25") at Side Deck
Colour(s)	Garnet Red (Matte) Euro Grey
Suspension	Enhanced Front and Rear Suspension
Heaviest Piece When Disassembled ⁴	18 kg
Mobility	
Turning Radius ⁴	1219 mm (48")
Max Speed ²	Up to 4.2 mph
Max. Obstacle Climb	75 mm (3")
Dual Breaking System	Regenerative and Electromechanical
Battery	
Battery Requirements ⁴	Lithium Ion Pack 8 Ah or 16 Ah MAXX Battery Pack
Battery Charger	3.5 - Amp Off-Board (lithiumion)
Total Weight Without Batteries ⁴	42kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or -) 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.

