

**GO GO**  
ENDURANCE **Li**



Colours Available:



# Travelling Has Never Been Easier!

## Enjoy a Smooth, Comfortable Ride

With the Go-Go® Endurance Li, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8 Ah airline-compliant lithium battery pack\*\*, the Go-Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need when you arrive at your destination. Enjoy a smooth, comfortable ride with Comfort-Trac Suspension.

Specifications	
<b>Scooter</b>	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length	1070 mm (42")
Overall Width	521mm (20.5")
Overall Height	883 - 920 mm (34.75" - 36.25")
Drive Tyres	Front: 50 x 203 mm (2" x 8") - Solid Drive Wheels: 67 x 203 mm (2.5" x 8") - Solid
Ground Clearance	83 mm (3.25") at Side Deck 36 mm (1.5") at Motor
Colour(s)	Garnet Red (Matte)
Suspension	Enhanced Front and Rear Suspension
Heaviest Piece When Disassembled	17.75 kg
<b>Mobility</b>	
Turning Radius	1225 mm (48.25")
Max Speed	4.2 mph (6.7 km)
Max. Obstacle Climb	40 mm (1.5")
Drive System	Rear-Wheel Drive 24v DC Motor
Dual Breaking System	Regenerative and Electromechanical
<b>Battery</b>	
Battery Requirements	Lithium Battery Box 8 Ah or 16 Ah
Battery Charger	24v 3.5A Off-Board (Lithium)
Total Weight Without Batteries	41 kg
Range Per Charge	8 Ah Battery: Up to 5 miles (8 km) at 23 st 3 lbs (147.5 kg) Up to 6.6 miles (10.6 km) at 14 st 4 lbs (91 kg) 16 Ah Battery: Up to 9.9 miles (16 km) at 23 st 3 lbs (147.5 kg) Up to 13.2 miles (21.3 km) at 14 st 4lbs (91 kg)

## Features



Dual Charging Ports Can Charge Battery Pack On/Off-board.



Delta Tiller With Ergonomic Wraparound Handles



Standard 8 Ah Airline-Compliant Lithium Battery