



Colours Available:



# Travelling Has Never Been Easier!

## Enjoy a Smooth, Comfortable Ride

With the Go Go Endurance Li, traveling has never been easier! This 4-wheel scooter has everything you need for a safe and convenient trip. With the standard 8 Ah airline-compliant lithium battery pack, the Go Go Endurance is perfect for air travel. Feather-touch disassembly provides the portability you need when you arrive at your destination. Enjoy a smooth, comfortable ride with Comfort-Trac Suspension.

### Specifications

Specifications	
Scooter	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length	1070 mm (42")
Overall Width	690 mm (27")
Overall Height	883 - 920 mm (34.75" - 36.25")
Drive Tyres	Front: 203 mm (8") - Solid Rear: 203mm (8") - Solid
Ground Clearance	83 mm (3.25") at Side Deck 36 mm (1.5") at Motor
Colour(s)	Garnet Red (Matte) Euro Grey (Matte)
Suspension	Enhanced Front and Rear Suspension
Heaviest Piece When Disassembled	17.75 kg
Mobility	
Turning Radius	1225 mm (48.25")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	40 mm
Drive System	Rear-Wheel Drive 24v DC Motor
Dual Breaking System	Regenerative and Electromechanical
Battery	
Battery Requirements	8 Ah (Standard) 16 Ah (Optional)
Battery Range	8 Ah 6.6 miles (10 km) 16 Ah 13.2 miles (21.3 km)
Battery Charger	24v 3.5A Off-Board (Lithium)
Total Weight Without Batteries	41 kg
Battery Weight	8 Ah: 3.5kg 16 Ah: 5.5kg

### Features



Dual Charging Ports Can Charge Battery Pack On/Off-board.



Delta Tiller With Ergonomic Wraparound Handles



Standard 8 Ah Li flight-safe battery



REV 02 - 30/01/2023

**PRIDE**  
MOBILITY PRODUCTS LTD.  
LIVE YOUR BEST