

FINESSE SPORT



Colours Available:



Viper Blue & Candy Apple Red

Sporty 8 mph Mobility Scooter

You're active and always on the go. You need a scooter that is up to the task! Appeal to your sporty side with the Apex Finesse Sport. Enjoy premium suspension, bright lighting, all in one sporty package.

Features



Cabin Light



Easy to Use Controls



Suspension

Specifications	
Scooter	
Weight Capacity	28 st 8lb (181 kg)
Overall Length ³	1190 mm (47")
Overall Width ³	648 mm (25.5")
Drive Tyres	254 mm
Ground Clearance ³	56 mm (2.2")
Suspension	Front & Rear
Heaviest Piece When Disassembled ³	35 kg
Colour(s)	Viper Blue & Candy Apple Red
Mobility	
Turning Radius ³	1566 mm (61.25")
Max Speed ¹	8 mph (12.8km/h)
Battery	
Battery Requirements ³	2 x 50 Ah Deep Cycle
Battery Charger	Off-Board
Battery Weight ⁴	14 kg Each
Total Weight Without Batteries ^{3,4}	73 kg

NOTE:

1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.

2 Includes standard seat weight.

3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.

4 Battery weight may vary based on battery manufacturer.

5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.