

GO GO
ELITE TRAVELLER *SPORT*



Colours Available:



Sleek, Sporty and Undeniably Stylish

The Go-Go® Elite Traveller Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter, as well as bright LED lighting.

Features



Feather-Touch Disassembly
With One Hand



Front and Rear LED Lighting



Extra 2" of Space For Added
Leg Room

Specifications	
Scooter	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length ³	1092 mm (43")
Overall Width ³	540 mm (21.25")
Drive Tyres	229 mm Front & Rear - Solid
Ground Clearance ³	38 mm (1.5") at Motor
Heaviest Piece When Disassembled ³	20 kg
Colour(s)	Red & Blue
Mobility	
Turning Radius ³	1283 mm (50.5")
Max Speed ¹	4 mph (6 km/h)
Battery	
Battery Requirements ³	2 x 18 Ah Deep Cycle 2 x 20 Ah Deep Cycle
Battery Weight ⁴	13 kg
Battery Charger	Off-Board
Total Weight Without Batteries ^{3, 4}	44 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3 %. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.