



Colours Available:




Revolutionary Suspension For The Smoothest Ride

There's smooth, and then there's the Apex Rapid with CTS smooth! This scooter has revolutionary suspension for the smoothest, most comfortable ride. From the park to your garden, this feature-rich scooter offers best-in-class ride quality.

Features




Compact, Foldable Seat
With Armrests



LED Lighting



Optional Colour Shrouds



Easy to Use Wrap Around
Delta Tiller

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length ³	1041 mm
Overall Width ³	521 mm
Drive Tyres	178 mm (7") Front - Solid 203 mm (8") Rear - Solid
Ground Clearance ³	38 mm
Heaviest Piece When Disassembled ³	16.6 kg
Colour(s)	Blue, Black, Red, Yellow, Orange & Pink
Suspension	Comfort-Trac Suspension - Front & Rear
Mobility	
Turning Radius ³	1162 mm
Max Speed ³	4 mph (6 km/h)
Battery	
Battery Requirements ³	2 x 12 Ah Deep Cycle
Battery Charger	2-Amp Off-Board
Total Weight Without Batteries ^{3, 5}	39.6 kg (87.4 lbs)
Battery Weight ⁴	18 Ah 13 kg

NOTE:

- 1 Varies with user weight, terrain type, battery amp-hour (AH), battery charge, battery condition, and tyre condition. This specification can be subject to a variance of + or - 10%.
- 2 Includes standard seat weight.
- 3 Due to manufacturing tolerances and continual product improvement, this specification can be subject to a variance of (+ or -) 3%.
- 4 Battery weight may vary based on battery manufacturer.
- 5 Does not meet airline requirements.

All specifications subject to change without notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice. Speed and range vary with user weight, terrain type, battery charge, battery condition and tire pressure. Due to manufacturing tolerances and continual product improvement, length, width, turning radius and ground clearance can be subject to variance of + or - 3%. Dimensions listed are for power base only. Overall measurements will vary based on seating and accessory selections.