

GO GO
ELITE TRAVELLER LX



Colours Available:



Sleek Style With Luxurious Features

Every day is a new adventure with the Go Go® Elite Traveller LX, an exceptionally stable and maneuverable travel scooter. With easy disassembly and a convenient storage basket, the Go Go Elite Traveller LX keeps up with your active life!

Specifications

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1060 mm (41.75")
Overall Width	521 mm (20.5")
Overall Height	885 mm (35")
Drive Tyres	Front: 178 mm (7") - Solid Rear: 203 mm (8") - Solid
Ground Clearance	36 mm (1.5") at Motor
Heaviest Piece When Disassembled	Front Section 17.75 kg
Suspension	Comfort-Trac Suspension
Colour(s)	Red & blue
Mobility	
Turning Radius	1219 mm (48")
Max Speed	4 mph (6.4km/h)
Drive System	Rear - Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Braking System	Regenerative and Eletromechanical
Maximum Obstacle Climbing Ability	50 mm
Battery	
Battery Requirements	2 x 12 Ah (Standard) 2 x 17 Ah (Optional)
Battery Weight	12 Ah - 9 kg 17 Ah - 14 kg
Battery Charger	2 - Amp Off-Board
Total Weight Without Batteries	41 kg

Features



Feather-Touch Disassembly
With One Hand



Large Back Pocket For Storage



Easy to Use Controls



REV 02 - 30/01/2023

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®