

Jazzy[®]
600 ES



Colours Available:



Extraordinary Manoeuvrability With a Flawless Design


The Jazzy[®] 600 ES has 14-inch drive wheels and independent front casters that handle terrain changes easily. Enjoy aggressive precision, a tight turning radius and a flawless design.

| Specifications | |
|--------------------------------|---|
| Chair | |
| Weight Capacity | 21 st 6 lbs (136 kg) |
| Overall Length | 1022 mm (40.2") |
| Overall Width | 622 mm (24.5") |
| Drive Tyres | 355 mm (14") Solid |
| Ground Clearance | 70 mm (2.75") at Battery Tray |
| Suspension | ATX Suspension (Active-Trac With Extra Stability) |
| Colour(s) | Red, Blue and Black |
| Power | |
| Turning Radius | 521 mm (20.5") |
| Max Speed | 4 mph (6 km/h) |
| Max. Obstacle Climb | 70 mm (2.75") |
| Battery | |
| Battery Requirements | 2 x 50 Ah, Deep Cycle |
| Battery Weight | 16.9 kg each |
| Battery Charger | 8 - Amp Off-Board |
| Total Weight Without Batteries | 79.1 kg (174.41lbs) |


Features



Fold Flat Seat



Extra Storage



Angle, Height and
Depth-Adjustable Armrest