



VR2 Controls



Lights & Indicators (6mph)



Solid or Pneumatic Tyre Options



Rear Suspension

# FUSION<sup>e</sup>

The Fusion E is a cost effective rear wheel drive power chair with a modular state of the art design that provides improved obstacle climbing performance and a smooth ride. The Fusion E is a versatile chair that will meet the demands of even the most active user.



## Specifications

### Power Chairs

Weight Capacity	23 stone 6 lbs (150 kg)
Total Weight without Batteries	86kg
Maximum Speed <sup>1</sup>	4 mph (6.4km/h)
Ground Clearance <sup>2</sup>	40 mm (1.57") at Lowest Point of Frame
Turning Radius <sup>2</sup>	641.5 mm (25.25")
Overall Size <sup>2</sup>	<b>Length:</b> 1050 mm (41.33") with Front Riggings <b>Width:</b> 630 mm (24.8")
Maximum Obstacle Climbing Ability	60 mm
Seat Sizes	15"W x 15"D to 22"W x 23"D
Controls	VR2 Controls
Angle Adjustable Backrest	- 5°, 15°, 25°, 35°

### Wheels

Rear Wheels	355 mm (14") Solid/Pneumatic/ Rear Wheel Drive
Front Wheels	254 mm (10") Pneumatic
Anti-Tip Wheels	50.8 mm (2")
Suspension	Independent Rear Suspension
Drivetrain	Two 350w Motor, Rear-Wheel Drive
Brakes	Auto Magnetic Brakes
Batteries <sup>3</sup>	2 x 50 Ah Deep Cycle
Battery Charger	5-Amp, Off-Board

### Other Features

Versatile Design, Easily Adjustable Seat to Allow for Growth

Swing-Away Legrests

Height Adjustable Armrests

- 1) Speed varies with user weight, terrain type, battery charge, battery condition and tire pressure.
- 2) Without front riggings.
- 3) Battery capacity may change when adding power positioning systems.
- 4) Weight capacity is dictated by the selection of seating.
- 5) All measurements +/- 1/8" with 18" x 18" seat, fixed towers, 21 stone 6 lbs. weight capacity and solid tires (add 3/8" with pneumatic tires). Widths may increase seat-to-floor height. Contact Customer Service for details. Listed measurements are in the lowest position; does not include elevating system unless specified.
- 6) Maximum dimensions of the Manual Recline vary by powerbase and seating system. Call for details.
- 7) Battery weight may vary +/- 3 lbs. based on manufacturer.
- 8) All Quantum Rehab® units can be ordered with specialty controls.
- 9) Calls for availability.

## Colours

Please contact your representative for the full range of options and colours.



**Black**