

ZOOF[®]



Zoof Classic Self Propelled

Colours Available:



Zoof Urban Transit

Colours Available:



Innovative and Adaptable Functionality in a Wheelchair

The Zoof mechanism makes storage and transportation a breeze, folding in one fluid motion with all parts integrated as one unit.

Specifications		
Scooter	Classic Self-Propelled	Urban Transit
Weight Capacity	21 st 2 lbs (135 kg)	21 st 2 lbs (135 kg)
Sitting Width	450 mm (17.7")	
Sitting Depth	480 mm (18.8")	
Sitting Height	530 mm (20.9")	
Overall Width	730 mm (28.7 ")	
Overall Length	1120 mm (44.1")	1020 mm (40.1")
Overall Height (push bar down)	980 mm (38.6")	
Width when Folded	600 mm (23.6")	
Length when Folded	480 mm (18.9")	440 mm (17.3")
Height when Folded	790 mm (31.1")	
Height of the Backrest	930 mm (36.6")	
Pivot Width	1150 mm (45.3")	
Turning Radius	1150 mm (43.3")	
Overall Weight	1 st 7 lbs (10.8 kg)	1 st 6 lbs (10.5 kg)
Wheel Sizes	Self Propel 24" Transit 12"	
Braking System	Stop Brake System	
Footrest	Mono frame	Double frame

Features



Stop Brake System



Quick release to switch between self-propel and transit wheels



Under Seat Storage



Foldable armrests to allow table use



Zoof Urban Transit
Flip up footplate



Zoof Classic Self-Propelled
Swing away footplate



REV 01 - 15/06/2023



Zoof Classic includes wheel bag with easy carry handles