

### The Zoof manual wheelchair (large wheels)

Model	Classic	Urban
Sitting width	450 mm / 17.7 in	
Sitting depth	480 mm / 18.8 in	
Sitting height	530 mm / 20.9 in	
Unfolded width	730 mm / 28.7 in	
Unfolded length	1120 mm / 44.1 in	1020 mm / 40.1 in
Unfolded height (push bar down)	980 mm / 38.6 in	
Folded width	600 mm / 23.6 in	
Folded length	480 mm / 18.9 in	440 mm / 17.3 in
Folded height	790 mm / 31.1 in	
Height of the backrest	930 mm / 36.6 in	
Pivot width	1150 mm / 45.3 in	
Weight without accessories	10,8 kg / 23.81 lbs	10,5 kg / 23.15 lbs
Weight large wheels	2,23 kg / 4.92 lbs each	
Weight brakes	0,28 kg / 0.62 lbs each	
Max. user weight	135 kg / 297.62 lbs	135 kg / 279.62 lbs
Max. weight on bag hook	5 kg / 11.02 lbs	
Max. weight on storage module	4 kg / 8.8 lbs	
Footrest	Mono frame	Double frame