

GO GO
ELITE TRAVELLER® 2

Our smallest turning
radius ever!



EZ TURN
TECHNOLOGY

Colours Available:



Scan the
QR code
for more
information!



Compact Turning Capabilities with extra style

The GoGo Elite Traveller 2 travels at speeds of up to 4 mph. With it's brand new turning radius technology at 37", it makes it our smallest turning radius yet. As well as its mesh chair making it one of the comfiest travel scooters on the market.

Specifications	
Scooter	
Weight Capacity	21 stone 6 lbs (136 kg)
Overall Length	1000 - 1036 mm (39.37 - 40.79")
Overall Width	494 mm (19.45")
Drive Tyres	Front: 178 mm (2" x 7") Solid Rear: 203 mm (2.5" x 8") Solid
Ground Clearance	34 mm (1.34")
Colour(s)	Tangerine orange, Black, Iceberg blue, Cardinal red, Robins egg blue, Sugar plum, Chell-O yellow
Mobility	
Turning Radius	940 mm (37")
Max Speed	4.23 mph (6.8 km/h)
Max. Obstacle Climb	50 mm (1.97")
Drivetrain	Rear-wheel drive, sealed transaxle, 24-volt DC motor
Dual Breaking System	Regenerative and electromechanical
Battery	
Battery Requirements	2 x 12-volt
Battery Range	4-wheel: (12 Ah) Up to 6.6 miles (10.62 km) (18 Ah) Up to 10.4 miles (16.7 km) Up to 13.8 miles (22.1 km) at 220 lbs. (99.79 kg)
Battery Weight	12Ah: 8.75 kg (19.5 lbs) 17Ah: 13.4 kg (29.6 lbs)
Total Weight Without Batteries	37.7 kg (83 lbs)

Features



37" turning radius (15% smaller than that of the Elite Traveller)



Easy 5 part Assembly



New stylish mesh seat for extra comfort



Delta tiller with high-level charge port



Arm rest storage bag comes as standard



REV 02 - 03/05/2024

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST[®]

www.pride-mobility.co.uk