

GOGO Mg²



Colours Available:



Ultra Lightweight foldable scooter

The GoGo Mg² travels at speeds of up to 4 mph. With its magnesium foldable frame at only 14.7 kg it makes it our lightest travel scooter yet and compact to fit in any car boot. With 10Ah lithium batteries as standard for a longer range per charge.

Specifications	
Scooter	
Weight Capacity	19 st 6lbs (124kg)
Overall Length	1016 mm (40")
Overall Width	429 mm (16.9")
Drive Tyres	Front: 203 mm (8") Solid Rear: 203 mm (8") Solid
Ground Clearance	74.9 mm (2.95")
Colour(s)	Black, Grey
Mobility	
Turning Diameter	2400 mm (94.5")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	8°
Drivetrain	Rear wheel
Dual Breaking System	Electromagnetic brake
Battery	
Battery Requirements	6Ah Li: (Optional) 10Ah Li: (Standard) 20Ah Li: (Optional)
Battery Weight	6Ah: 1 kg 10Ah: 1.7 kg 20Ah: 3 kg
Total Weight Without Batteries	14.7 kg

Features



Optional storage available



Lightweight lithium batteries.



Foldable compact design.



Optional arm rests available



REV 02 - 25/03/2024

