

Jazzy Ultra-Lite
Carbon

Weighing
only 14.5kg



Colours Available:



The Ultimate Lightweight Chair

Travel made easy

This carbon fibre power chair has plenty to offer with its one touch fold release it's never been easier to get moving, whilst weighing only 14.5 kg making it our lightest powered wheelchair yet.

Specifications	
Chair	
Weight Capacity	21 st 5 lb (136 kg)
Opened Dimensions	1000 x 620 x 950 (mm)
Folded Dimensions	620 x 340 x 800 (mm)
Drive Tyres	Front 203 mm (8") Rear 304 mm (12")
Colour(s)	Black
Power	
Turning Radius	900 mm (35.4")
Seat Width	450 mm (17.71")
Seat Depth	450 mm (17.71")
Max Speed	4 mph (6 km/h)
Battery	
Battery Range	10 Ah: 9 miles (14.8 km) 20 Ah: 18 miles (28.96 km)
Battery Requirements	10 Ah Li (standard) 20 Ah Li (optional)
Battery Weight	10 Ah: 1.4 kg (standard) 20 Ah: 3 kg (optional)
Battery Charger	29.4v 2A
Total Weight Without Batteries	14.5 kg


Features




Flight safe batteries



Fold away footplate



Compact one lever
fold function



Stretch mesh backrest for
lightweight extra comfort



REV 03 - 07/11/2025

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®