

GO•GO[®]
SUPER PORTABLE



Colours Available:



Super Lightweight folding scooter

The GoGo Super Portable travels at speeds of up to 4 mph. Weighing only 18.9kg its carbon fibre frame makes it a strong yet lightweight. As well as it's unique folding mechanism for easy transportation.

Super lightweight foldable scooter


GO GO
SUPER PORTABLE



Specifications	
Scooter	
Weight Capacity	21 st 6lbs (136kg)
Overall Length	945 mm (37.2")
Overall Width	547 mm (21.54")
Drive Tyres	Front: 178 mm (7") Solid Rear: 203 mm (8") Solid
Ground Clearance	47 mm (1.85")
Colour(s)	Black, Mint
Mobility	
Turning Radius	1311 mm (51.62")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	40mm
Dual Breaking System	Regenerative and Electromagnetic brake
Battery	
Battery Requirements	11.7 Ah Li: (Standard) Flight compatible 16 Ah Li: (Optional)
Battery Range	11.7 Ah Li: Up to 10.1 miles (16.3 km)
Battery Weight	11.7 Ah: 3.6 kg 16 Ah: 4.4 kg
Total Weight Without Batteries	18.9 kg




360 degrees magnetic charging point



Rear footpeg for upright storage option



Underseat Storage



Optional arm rests available



Front and rear lighting



USB device charging port.



Storage bag as standard



Adjustable front tiller