

GO GO
ELITE TRAVELLER *SPORT*



Colours Available:



Patented EZ Technology
offering 37" turning radius

Sleek, Sporty and Undeniably Stylish

The Go Go Elite Traveller® Sport offers simple disassembly and an ergonomic delta tiller, perfect for one-handed operation. Enjoy the versatility of this sporty scooter, as well as bright LED lighting.

Specifications	
Scooter	
Weight Capacity	23 st 3 lbs (147 kg)
Overall Length	1092 mm (43")
Overall Width	540 mm (21.25")
Overall Height	889 - 914 mm (35 - 36")
Drive Tyres	Front: 229 mm (9") - Solid Rear - 229 mm (9") - Solid
Ground Clearance	76.25 mm (3") at Mid-Frame 38 mm (1.5") at Motor
Heaviest Piece When Disassembled	Front Section: 20 kg
Colour(s)	Red & Blue
Mobility	
Turning Radius	1283 mm (50.5")
Maximum Obstacle Climbing Ability	50 mm (2")
Max Speed	4 mph (6.4km/h)
Drive System	Rear-Wheel Drive Sealed Transaxle 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Battery	
Battery Requirements	2 x 20 Ah
Battery Range	20 Ah 14.5 miles (23 km)
Battery Weight	12.5 kg
Battery Charger	20 Ah 2 A Off-Board Charger 20 Ah Extended Range 3.5A Off-Board Charger
Total Weight Without Batteries	44 kg

Features



Feather-Touch Disassembly
With One Hand



Front and Rear LED Lighting



Extra 2" of Space For Added
Leg Room



Patented EZ Technology
offering 52" turning radius