



# Ultra Lightweight portable scooter

The GoGo Mg lite travels at speeds of up to 4 mph. With its magnesium foldable frame at only 12 kg and compact design it fits in any car boot. All with 11Ah lithium batteries making it your flightsafe travel companion.

Specifications	
<b>Scooter</b>	
Weight Capacity	18 st 9lbs (120kg)
Overall Length	990 mm (38.9")
Overall Width	625 mm (24.6")
Drive Tyres	Front: 195 mm (7.5") Solid Rear: 195 mm (7.5") Solid
Ground Clearance	74.9 mm (2.95")
Colour(s)	Black
<b>Mobility</b>	
Turning Radius	1450 mm (57")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	25 mm
<b>Battery</b>	
Battery Requirements	11 Ah Lithium (flightsafe)
Battery Range	9.3 miles (15 km)
Battery Weight	1.6 kg
Total Weight Without Batteries	12 kg

## Features



Digital Display with contactless lock



Battery Level Display.



Easy disassembly



Lightweight Quick Release Seat with adjustable armrests



Front and rear brake light



REV 02 - 03/07/2025

**PRIDE**  
MOBILITY PRODUCTS LTD.  
LIVE YOUR BEST®

[www.pride-mobility.co.uk](http://www.pride-mobility.co.uk)