

GO CHAIR



Colours Available:



Enhance Your Independence! And Enjoy Your Favourite Places

The Go Chair® is number one for portability! With feather-touch disassembly, you can take your Go Chair to all your favourite places and enjoy greater independence.

Specifications	
Chair	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	842 mm (33")
Overall Width	575 mm (22.5")
Front/Rear Tyres	Front: 152mm (6") - Solid Rear: 229 mm (9") - Solid
Ground Clearance	45mm (1.75")
Heaviest Piece When Disassembled	16.5 kg
Colour(s)	Red, Orange, White, Blue, Yellow, Pink
Power	
Turning Radius	648 mm (28.5")
Max Speed	Up to 4 mph (6.4 km/h)
Drivetrain	Rear-Wheel Drive, Two Motor
Dual Breaking System	"Intelligent Braking" Electronic Regenerative, Disc Park Brake
Battery	
Battery Requirements	2 x 18 Ah Batteries
Battery Range	Up to 13.2 miles (21.24 km)
Battery Weight	14 kg
Battery Charger	2 - Amp, Off-Board
Total Weight Without Batteries	44.25 kg

Features



Disassembles into Five Pieces



Large Comfortable Foot Plate



Adjustable Height, Width & Armrest



REV 03 - 03/06/2024

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®