



Colours Available:



Ultra Lightweight travel scooter

Enjoy worry-free travel with the GoGo Ultra, a compact scooter that is simple to disassemble. Ride in comfort additional leg room and new delta tiller.

Specifications	
Scooter	
Weight Capacity	18 st 9lbs (120kg)
Overall Length	1022 mm (40.23")
Overall Width	490 mm (19.29")
Drive Tyres	Front: 203 mm (8") Solid Rear: 203 mm (8") Solid
Ground Clearance	32 mm (1.25")
Colour(s)	Matte Black, Matte Grey
Mobility	
Turning Radius	1385 mm (54.5")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	57.5mm (2.25")
Battery	
Battery Requirements	2 x 12 ah Sealed Lead Acid
Battery Range	11.5 miles (18.5 km)
Battery Weight	4 kg
Total Weight Without Batteries	36.3 kg

Features



New longer footplate for additional leg room



New delta tiller



8" alloy wheels



Front storage basket



REV 01 - 30/06/2025

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®