

Jazzy Carbon

SWIFT



Weighing only
10.26 kg!

Extra lightweight innovative powerchair

The Jazzy Carbon Swift is a lightweight powerchair weighing just 10.26 kg, making it perfect for travel and everyday use. With easy-to-use controls, speeds of up to 4 mph, and a range of up to 9.3 miles, it offers reliable mobility in a compact, portable design.

Extra lightweight innovative powerchair

Specifications	
Scooter	
Weight Capacity	18 st 10lbs (120kg)
Overall Length	800 mm (31.4")
Overall Width	565 mm (22.24")
Drive Tyres	Front: 125 mm (5") Solid Rear: 165 mm (6.5") Solid
Ground Clearance	74.9 mm (2.95")
Colour(s)	Black
Mobility	
Turning Radius	900 mm (35.4")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	25 mm
Battery	
Battery Requirements	10 Ah Lithium (flightsafe)
Battery Range	9.3 miles (15 km)
Battery Weight	1.6 kg
Total Weight W/O batteries and cushion	10.26 kg


Features



Narrow frame offering effective maneuverability



Flight safe batteries



Easy-fold mechanism



Easy-to-use controls



REV 01 - 30/06/2025

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®

www.pride-mobility.co.uk