



Weight capacity of
up to 28.7st!



Heavyduty Powerchair lightweight frame

The Jazzy Carbon Ultra Max, a revolutionary heavy-duty powerchair that combines exceptional strength with ultra-lightweight design. Engineered for performance and comfort, it supports users up to 28.7 stone while maintaining an incredibly light frame of just 11.83kg.

Specifications	
Scooter	
Weight Capacity	28 st 7lbs (182 kg)
Overall Length	990 mm (38.9")
Overall Width	625 mm (24.6")
Drive Tyres	Front: 175 mm (7") Solid Rear: 210 mm (8) Solid
Ground Clearance	75 mm (3")
Colour(s)	Blue
Mobility	
Turning Radius	850 mm (25.5")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	25 mm
Battery	
Battery Requirements	12 Ah Lithium (flightsafe)
Battery Range	9.3 miles (15 km)
Battery Weight	2.1 kg
Total Weight W/O batteries and cushion	11.83 kg

Features



Heavy duty weight capacity of up to 28.7 st



Lightweight strong carbon fibre frame



Folding footplate



Easy folding mechanism



REV 01 - 30/06/2025

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®

www.pride-mobility.co.uk