



# Lightweight Aluminium travel scooter

Experience the convenience of a lightweight, compact travel scooter with airline safe batteries. With its new Delta Tiller and additional legroom the Aluma-Lite 2.0 is our most advanced aluminium scooter yet!

Specifications	
<b>Scooter</b>	
Weight Capacity	18 st 9lbs (120kg)
Overall Length	1022 mm (40.2")
Overall Width	490 mm (19.29")
Drive Tyres	Front: 203 mm (8") Solid Rear: 203 mm (8") Solid
Ground Clearance	32 mm (1.25")
<b>Colour(s)</b>	Matte black, Matte grey
<b>Mobility</b>	
Turning Radius	1385 mm (54.5")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	57.5 mm (2.25")
<b>Battery</b>	
Battery Requirements	12 Ah Lithium (flightsafe)
Battery Range	11.5 miles (18.5 km)
Battery Weight	2.5 kg
Total Weight Without Batteries	31.8 kg

### Features



New delta tiller



New extended footplate for extra leg room.



New mesh seat for added comfort



8" alloy wheels



REV 01 - 30/06/2025

