



Colours Available:



Revolutionary Suspension For The Smoothest Ride

There's smooth, and then there's the Apex Rapid with CTS smooth! This scooter has revolutionary suspension for the smoothest, most comfortable ride. From the park to your garden, this feature-rich scooter offers best-in-class ride quality.

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1060 mm (41.75")
Overall Width	521 mm (20.5")
Drive Tyres	Front - 203 mm (8") - Solid Rear - 203 mm (8") - Solid
Ground Clearance	36 mm (1.5") At Motor
Heaviest Piece When Disassembled	Front Section: 17.75 kg
Colour(s)	Yellow, Blue, Black, Orange, Pink, Red
Suspension	Comfort-Trac Suspension - Front & Rear
Mobility	
Turning Radius	1219 mm (48")
Max Speed	4 mph (6.4 km/h)
Drive System	Rear - Wheel Drive, Sealed Transaxle, 24v DC Motor
Dual Braking System	Regenerative and Electromechanical
Maximum Obstacle Climbing Ability	50 mm
Battery	
Battery Requirements	2 x 12 Ah (Standard) 2 x 18 Ah (Optional)
Battery Range	12 Ah: 5.5 miles (8.75 km) 18 Ah: 10.5 miles (17 km)
Battery Charger	2 A Off-Board Charger
Total Weight Without Batteries	41 kg
Battery Weight	12 Ah: 9kg 17 Ah: 14kg

Features




Compact, Foldable Seat With Armrests



LED Lighting



Optional Colour Shrouds



Easy to Use Wrap Around Delta Tiller