

ATMOS



Powerful Performance and built to last.

Enjoy premium suspension, bright lighting, a weight capacity of up to 25st, speeds of up to 8mph and other sophisticated features, all in one robust package.

Specifications

Specifications	
Scooter	
Weight Capacity	25st (159 kg)
Overall Length	1420 mm (55.9")
Overall Width	687-707 mm (27-27.83")
Drive Tyres	Front: 330 mm (13") - Pneumatic Rear: 330 mm (13") - Pneumatic
Heaviest Piece When Disassembled	91.26 kg
Colour	Blue
Suspension	Front & Rear Suspension
Mobility	
Turning Radius	2090 mm (82.28")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	100 mm
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 75 Ah Deep Cycle
Battery Range	Up to 27 miles (43.5 miles)
Battery Charger	Off-Board
Battery Weight	22.14 kg
Total Weight Without Batteries	113.39 kg


Features



Front & Rear Suspension



Generously sized footplate for extra room



Adjustable armrests