



Colours Available:



Bold, Adventurous and Innovative

The BAJA Bandit is the first of the BAJA range to hit the UK market and has a whole new adventure in mind. The Bandit is great for rough terrain, with great suspension and an innovative breaking system. Travelling off-road at speeds of up to 8mph has never been easier! Furthermore with lithium batteries this BAJA powerhouse will keep going all day long at no extra weight.

Specifications	
Scooter	
Weight Capacity	28 st 8 lbs (181 kg)
Overall Length	1260 mm (49.6")
Overall Width	707 mm (27.83")
Front/Rear Tyres	Front: 239 mm (9.4") Solid Rear: 239 mm (9.4") Solid
Ground Clearance	40 mm (1.57")
Colour(s)	Tan, Sage
Mobility	
Turning Radius	1347 mm (53")
Max Speed ¹	8 mph (12.9 kmh)
Battery	
Battery Requirements	20 Ah Li 40 Ah (optional) Li
Battery Range	20 Ah: Up to 14.5 miles (23.34 km) 40 Ah: Up to 29 miles (46.72 km)
Battery Weight	6.62 kg each
Battery Charger	4-amp Off-Board
Total Weight Without Batteries	72.3kg

Features



Added suspension for a smooth ride even over rough terrain



Lithium batteries for longer range



Compact frame for easy manoeuvrability and turning radius of only 53"



9.4" Front and rear tyres