



Colours Available:



Style and Comfort At 6 mph

Embrace your potential! The Colt Deluxe 2.0 is fully loaded with style and advanced features, including frosted LED headlights and an integrated storage space with cup holder. Drive at speeds up to 6 mph and see what life has in store for you!



| Specifications | |
|----------------------------------|---|
| Scooter | |
| Weight Capacity | 28 st 8 lbs (181 kg) |
| Overall Length | 1194 mm (47") |
| Overall Width | 560 mm (22") |
| Drive Tyres | 254 mm Front & Rear (10")- Solid |
| Ground Clearance | 64 mm (2.5") |
| Heaviest Piece When Disassembled | Front Section 29.5 kg |
| Colour(s) | Red, Matte Blue, Grey |
| Suspension | Front & Rear Suspension Front & Rear Suspension |
| Mobility | |
| Turning Radius | 1334 mm (52.5") |
| Max Speed | 6 mph (10 km/h) |
| Max. Obstacle Climb | 80 mm (3.25") |
| Drivetrain | 24v DC Motor |
| Dual Breaking System | Regenerative & Electromechanical |
| Battery | |
| Battery Requirements | 2 x 40 Ah Deep Cycle |
| Battery Range | Up to 21 miles (39.2 km) |
| Battery Charger | 5-A Off-Board |
| Battery Weight† | 14.75 kg |
| Total Weight Without Batteries | 68 kg |

Features



Comfortable High Back Seat



Integrated USB Charging Point for Device



Pneumatic Tyres With Solid Tyres Available



REV 04 - 31/05/2024

