



Colours Available:



# Robust Features With Supreme Quality

The Colt Pursuit drives at speeds of up to 8 mph. As well as being equipped with a regenerative braking system, a comfortable high-back seat and complete LED lighting



Specifications	
<b>Scooter</b>	
Weight Capacity	25 st (159 kg)
Overall Length	1353 mm (53")
Overall Width	641 mm (25")
Drive Tyres	Front: 254 mm (13") Pneumatic Rear: 254 mm (13") Pneumatic
Ground Clearance	89 mm (3.5")
Heaviest Piece When Disassembled	19 kg (seat)
Colour(s)	Silver & Black
Suspension	Front & Rear Suspension
<b>Mobility</b>	
Turning Radius	1695 mm (67")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	100 mm
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
<b>Battery</b>	
Battery Requirements	2 x 50 Ah (standard) 2 x 70 Ah (optional)
Battery Range	Up to 34.8 miles (56 km)
Battery Weight	50Ah : 14.8 kg 70Ah: 24 kg
Battery Charger	8-Amp Off-Board
Total Weight Without Batteries	78.5 kg

### Features



High Visibility Front and Rear LED Lighting



Indicators



13" Alloy Style Pneumatic Tyres



REV 03 - 31/05/2024

