



Colours Available:



Sporty Style with Advanced Comfort

The all new Colt Sport is fully loaded with sporty style and advanced features, including a new DX Tiller, LED headlights and an advanced suspension seat post as well as an all new memory foam adjustable seat. Experience all of this whilst travelling speeds of up to 8 mph.

Specifications	
Scooter	
Weight Capacity	25 st (159 kg)
Overall Length	1194 mm (47")
Overall Width	657 mm (26")
Drive Tyres	Front: 254 mm (10") Solid or Pneumatic Rear: 254 mm (10") Solid or Pneumatic
Ground Clearance	57 mm (2.25")
Heaviest Piece When Disassembled	29.75 kg (65.5lbs)
Colour(s)	Matte Black and Stone Grey
Suspension	Front & Rear Suspension
Mobility	
Turning Radius	1334 mm (52.5")
Max Speed	8 mph (12.9 km/h)
Max. Obstacle Climb	50 mm (2")
Drivetrain	24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Battery	
Battery Requirements	2 x 40Ah (standard) 2 x 50Ah (optional)
Battery Range	Up to 16.7 miles (26.8 km)
Battery Weight	40 Ah: 14.75kg 50 Ah: 16.5 kg
Battery Charger	Off-Board
Total Weight Without Batteries	68 kg

Features



Suspension Seat Post



10" Alloy Style Pneumatic or solid Tyres (solid tyres only available in black)



Memory foam adjustable seat



DX Tiller and LED Lighting