

MAXIMA



Colours Available:




Ultra Heavy Duty With Maximum Performance

Heavy duty and robust, the Maxima scooter maximises your potential! Drive at speeds up to 6 mph and experience a fully adjustable tiller. Enjoy the convenience of the front basket, perfect for storing your personal items.

(Specifications		
Scooter	4-Wheel	3-Wheel
Weight Capacity	35 st 10 lbs (226.75 kg)	
Overall Length	1206 mm (47.5")	
Overall Width	686 mm	
Drive Tyres	Front: 265 mm (10.4") - solid Rear: 270mm (10.75") - solid	Front: 257 mm (10.1") - solid Rear: 270mm (10.75") - solid
Ground Clearance	63.5 mm (2.5")	57.15 mm (2.25")
Colour(s)	Blue (Matte), Red (Matte)	
Mobility		
Turning Radius	1275 mm (49")	1346 mm (53")
Max Speed	6 mph (9.6 km/h)	
Seat-to-Ground Height With Power Elevating Seat	561.34 mm - 612.14 mm	58.93 mm - 725.68 mm with Power Elevating Seat
Seat-to-Deck Height With Power Elevating Seat	435 mm - 485.9 mm	440. 94 mm - 600 mm
Drivetrain	24v Single Motor	
Dual Breaking System	Regenerative & Electromechanical	
Battery		
Battery Requirements	2 x 50 Ah	
Battery Range	Up to 19.8 miles (13.8 km)	Up to 19.5 miles (31.4 km)
Battery Charger	5A Off-Board	5A Off-Board
Total Weight Without Batteries	83 kg	89 kg
Battery Weight	17 kg	


Features



Easy Control Tiller



Available in a 3 or 4 Wheeler

Comfortable Medium
- Back Seat

REV 04 - 03/06/2024

PRIDE
MOBILITY PRODUCTS LTD.
LIVE YOUR BEST®

www.pride-mobility.co.uk