



Colours Available:



Lightweight Aluminium Travel Scooter

One of the lightest aluminum travel scooter on the market! Experience the convenience of a lightweight, compact travel scooter with airline safe batteries. Easy to put into the boot, with its heaviest piece weighing only 11.6 kg. The Apex Aluma-Lite plus, with its aluminum frame, weighs only 33.5kg making it one of our lightest four wheeled scooter yet!

Specifications	
Scooter	
Weight Capacity	21 st 6 lbs (136 kg)
Overall Length	1010 mm (39.76")
Overall Width	495 mm (19.48")
Overall Height	840 mm (33")
Drive Tyres	Front: 203 mm (8") Solid Rear: 203 mm (8") Solid
Ground Clearance	32 mm (1.25") at Motor
Heaviest Piece When Disassembled	11.6 kg - Rear Main Frame
Colour(s)	Red, Blue & Black
Mobility	
Turning Radius	1385 mm (54.5")
Max Speed	4 mph (6.4mph)
Max. Obstacle Climb	57.5 mm (2.25")
Drivetrain	Rear-Wheel Drive Single Motor 24v DC Motor
Dual Breaking System	Regenerative & Electromechanical
Maximum Obstacle	57.5 mm (2.25")
Battery	
Battery Requirements	12 Ah (Standard) 30 Ah (Optional)
Battery Range	12 Ah: 8.9 miles (14.4 km) 30 Ah: (Optional): 25.5 miles (41 km)
Battery Charger	Off-Board
Total Weight Without Batteries	33.5 kg
Battery Weight	12 Ah: 4.25 kg 30 Ah 8.63 kg

*30 Ah Battery is not a flight safe battery

Features



Suspension



Easy to Use Controls



REV 03 - 31/05/2024

