



Colours Available:



Our Lightest Aluminium Travel Scooter

One of the lightest aluminum travel scooter on the market! Experience the convenience of a lightweight, compact travel scooter with airline safe batteries. Easy to put into the boot, with its heaviest piece weighing only 11.2 kg. The Apex Aluma-Lite, with its aluminum frame, weighs only 28.8 kg making it our lightest four wheeled scooter yet!

Specifications	
Scooter	
Weight Capacity	18 st 13 lbs (120 kg)
Overall Length	1022 mm (40.25")
Overall Width	490 mm (19.25")
Front/Rear Tyres	Front: 203 mm (8") Solid Rear: 203 mm (8") Solid
Ground Clearance	32 mm (1.25")
Heaviest piece when disabled	11.2 kg -
Colour(s)	Red, Blue & Black
Mobility	
Turning Radius	1385 mm (54.5")
Max Speed ¹	4 mph (6.4 kmh)
Battery	
Battery Requirements	12 Ah Lithium Battery
Battery Range	8.9 miles (14.4 km)
Battery Weight	2.5 kg
Battery Charger	2-amp Off-Board
Total Weight Without Batteries	28.8kg

Features



Easy-to-Grip Tiller
Adjustment Knob



Airline Safe Batteries



Lightweight
Aluminum Frame



8" Front Rear Wheel