

Jazzy

CARBON MAX HD

Extra wide
20" seat



Heavyduty Powerchair lightweight frame

The Jazzy Carbon Ultra Max HD, a revolutionary heavy-duty powerchair with additional seat width of 19.5". With all the unique features of the Jazzy Carbon HD. Exceptional strength with ultra-lightweight design. Engineered for performance and comfort, it supports users up to 28.7 stone while maintaining an incredibly light frame of just 13.8kg.

Specifications	
Scooter	
Weight Capacity	28 st 7lbs (182 kg)
Overall Length	950 mm (37.4")
Overall Width	665 mm (26.8")
Drive Tyres	Front: 175 mm (7") Solid Rear: 270 mm (10.6) Solid
Ground Clearance	75 mm (3")
Colour(s)	Red
Mobility	
Turning Radius	925 mm (36.4")
Seat Width	510 mm (20")
Seat Depth	470 mm (18.5")
Max Speed	4 mph (6.4 km/h)
Max. Obstacle Climb	25 mm
Battery	
Battery Requirements	20 Ah Lithium
Battery Range	12.6 miles (20 km)
Battery Weight	3.4 kg
Total Weight W/O batteries and cushion	13.8 kg

Features



Heavy duty weight capacity of up to 28.7 st



Heavy duty seat 19.5" wide



Folding footplate



Easy folding mechanism



REV 01 - 17/02/2026

